

## Cleaning and maintenance

- You should take out the plug after you use, let it cool completely then start to clean it.
- Please use wet cloth clean the dirt on the product after you use it.
- Please don't use the corrosive cleaning fluid, benzene or diluents to clean the product.
- Avoid to trouble, please don't use wet cloth to clean the control panel, you can only use dry cloth to clean it.
- If you don't use it for a long time, please clean it and put in carton, and put in a ventilated and dry place, in case the product get damp and affect the performance of the product.

## Trouble shooting

This chapter summarizes the most common problems encountered in the use of this product. If you can't solve the problem with the following information, please consult with the maintenance center or customer service center authorized by your district.

If you encounter the other questions didn't mention as the list, please try as below steps:

1. Turn off the machine, take out the power plug.
2. Keep the products in off state 60 minutes.
3. Plug the power plug into the power socket again and start the machine.

Problem	Possible causes	Solution
The product can't work	The power plug is out or adaptor DC end is not plug in.	Plug the power plug in the power socket, or insert the adaptor DC end to the product, and check the product voltage is it well matched with local voltage.
	After 10 minutes massage, the product switched to standby mode automatically.	If you want to massage again, press the "on/off" button to start.
	You continuous use the product for 40 minutes, and the overheat protection device cut off the power.	You must take out the plug and cool for 60 minutes if you want to reuse it.
The product get heat	Overload on the product (Somebody stand on the products or with weights on the product)	Take out the power plug, move out the weights (if the weights still on the products), plug in the power plug again and start the machine.
	The warm function is on	Turn off the warm function
With noise come out	You continuous use the product for 40 minutes.	Take out the power plug and let it cool, take out the power plug every time after you use it.
	The problem caused by the mechanical structure.	It is normal noise, not mean the product has function problem.

## Product operation sketches



Kneading massage on neck, shoulder and back

## Product specification

Model: SL-D501

Name: Massage Belt

Rated voltage: 12V

Rated power: 25W


Rated working time: 10 minutes

Safety structure:

**iRest**<sup>®</sup>

## U SHOULDER HUG Stretch Your Neck And Shoulder

SL-D501



**SAVE THESE INSTRUCTIONS**

Operation Instructions

※Products refer to the actual product, without prior notice.

## Safety And Maintenance



### 1. Important Safety Warning

- Do not allow any foreign matter to enter the clearance between the massage head and the casing
- Do not allow the product to come into contact with water so as to avoid any possible accidents that may result from a short circuit.
- If the mains cable or plug is damaged, do not disassemble and attempt to repair it by yourself. Please contact your distributor.
- Ensure that the mains power is suitable for this product.
- Please disconnect the product from the mains supply after use or before cleaning to avoid damage to the equipment or personal injury.
- Please use this product according to this Manual.
- Do not use any spare parts and attachments other than those recommended.
- Do not use this product if its casing or mains cable is damaged.
- Do not use this product if you are under the influence of alcohol or feeling unwell.
- Do not use this product with too much force so as to avoid injury.
- Do not keep this product on for more than 15 minutes so as to avoid overheating.
- Do not use this product within one hour of eating.
- Please use this product carefully according to your doctor's advice if you have a neck or back injury or if you have had an operation recently.

### 2. Working Environment

- Do not use this product in a highly humid environment such as a bathroom.
- Do not use this product immediately after there has been a sharp change in the environmental temperature.
- Do not use this product in an environment where there are lots of dust or corrosive gases.
- Do not use this product when it is covered by a blanket or pillow, etc.

### 3. People That Should Not Use This Product

- People with contagious skin conditions
- Patients with osteoporosis
- Patients with heart diseases and with electronic medical instruments like a heart pacemaker
- People in a high fever
- Women that are pregnant or menstruating
- Those who are injured or suffering from any body surface disease
- The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.
- Children under 14-year old or people mentally unbalanced without being Supervised are not allowed to use this machine.
- Those who to need rest according to a doctor's advice or those who do not feel well
- Do not use this product when you are wet
- Do not use this product on a pet.

### 4. Safety Precautions


- Do not plug or unplug this product when your hands are wet
- Do not pull the cable when unplugging. Do not plug or unplug carelessly.

## Safety And Maintenance



- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- Do not damage the cable or change the circuit of the product.
- Please cease using this product in case of a power cut in order to avoid any body injury upon sudden power resumption.
- Please stop using this product and consult your local supplier if you notice anything abnormal with the product during use.
- Please stop using this product and consult your doctor for advice
- If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
- The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance.
- Forbidden to use when the cover broken.
- Forbidden to press it with stronger force or by whole body, don't touch it with hard objects, don't put it under heavy stress or tread on it.
- Prevent water or other corrosive liquid inflow to the products and cause stoppage and damage the machine.
- Please don't touch the wire line or controller line with hand, it may be cause the machine suspension and swing, lead to stop working because of the wire core was broken.
- Please don't change the structure of the machine or replace the spare parts.
- Keep it away from fire or inflammable place, it is more prone to stoppage and danger, cause distortion, discolored or damaged.
- If any abnormal happens when using, please stop it and cut off the power.
- People who mentally unbalanced or paralyzed should be carefully to use and under the watchful eye of others or consult to the doctors before use.


## Name and function of components



## Function

- Fashionable and streamline appearance patent design, novelty and with science and technology feeling.
- Strong structure and fashionable color design. Stronger toughness and better in resistance.
- The kneading head covered with antibacterial and pro-skin fabric, smooth surface and good resilience, make the massage experience softer and more comfortable.
- Unique 3D humanoid massage design, according with human body engineering, adjust the massage compactness freely, easy to operate.
- 2 speeds and 2 modes freely for choice, rated working time is 10 minutes, stop automatically, safety and convenient to use.
- With high stable DC motor, run quieter.
- With overload current protection function, anti-interference and anti-static function.
- Promote blood circulation on neck, shoulder and back area.

## Method of usage



- Insert the power line plug to the product connector. (Make sure they are connected well, avoid to falling out)
- Insert the adaptor into the power socket. When the red indicator light twinkles, means the power on. (Make sure they are connected well, avoid to falling out)
- Press the "on/off" button, the red indicator light on, start the manual mode automatically, "H" indicator light on, with "H"(high speed) to massage
- Press the "manual" button, you can switch high speed and low speed, and the corresponding button light "H" and "L" will be lightened.
- Press the "auto" button to start auto mode, and the "I" indicator light will be lightened. This button is switchable between "I" and "II", and the corresponding button light "I" and "II" will be lightened. The product will stop to running after 10 minutes.
- You can choose warm function when you massage. Press "warm" button, the indicator light will be lightened, the massage head part start to radiates heat. Press this button to on/off warm function.
- Press "on/off" button to finish all of the function.